

IndiRay®

Medical X-Ray Film Viewer - LED

True 10,000 LUX for CT, MRI



Uniform Clear Vision

No drop in Light Intensity

IndiRay® Medical X-Ray Film Viewer - LED



Unique Features

- No drop in light intensity. Dissipation of heat technology enhances life of LED.
- Critical components from Japan and Germany ensure reliable performance
- Film insert activated switch for each panel
- Continuous brightness control for each panel
- Low power consumption
- Sleek design. Thickness: 25 mm only. Brush finished anodized Aluminum frame. Provides elegant photo frame look.
- Lightest in weight
- Mounting mechanism through covered holes provided at four corners. Restricts dust flow, Ideal for OT

Model	Viewing Area cm	Luminous Density approx. Lux*	Light Source & Temp	Power Consumption Watts	Power Input	Outer Dim. L x H x B (cm)
Single Film HLX 1	44 x 36	10000	LED 6500 K	15	220-240 VAC	44 x 52 x 2.5
Two Films HLX 2	44 x 72	10000	LED 6500 K	30	220-240 VAC	82 x 52 x 2.5
Three Films HLX 3	44 x 108	10000	LED 6500 K	45	220-240 VAC	120 x 52 x 2.5
Four Films HLX 4	44 x 144	10000	LED 6500 K	60	220-240 VAC	160 x 52 x 2.5

Product specifications and features are subject to change



Sterling Imaging Solutions
Mumbai, India
E: sterling@sterlingimaging.com | W: www.sterlingimaging.com