



Vertical Bucky Stand

Elegant and compact design in both floor and wall mounted models, which meet the needs of every user for upright radiography. Fast and simple routine x-ray exposures for a standing and a sitting patient.

Features:

- Perfect counterbalancing for smooth vertical movement with least efforts.
- Simple, smooth and effortlessly adjustable to any desired height.
- Integrated chin recess that provides comfort to the patient and convenient patient positioning.
- Simple and quick installation – ready to operate in any radiological department.
- Bucky loading from both sides – left or right.

Specifications:

Standard type	Non-tilting bucky
Bucky	Par speed
Max. cassette size accommodation	430 mm x 430 mm or 17" X 17"
Vertical bucky travel	1370 mm (from 380 to 1750 mm)
Film centre to floor min. distance	380 mm
Front top to film distance	approx. 40 mm
Vertical movement	Manual, fully counterbalanced
Locks	Mechanical breaks (Electromagnetic optional)
Electrical requirement	Single phase 230V +/- 10%. 50/60 Hz
Dimension	620 W x 625 D x 2070 H (mm)
Weight	75 Kg. approx.
Grid (Standard)	17 1/4 x 18 7/8, 85 lines, 8:1 ratio, FD 150 cm