


IndiRay®

Medical X-Ray Film Viewer - LED

CT MRI X-Ray



Uniform Clear Vision

No drop in Light Intensity

IndiRay[®] Medical X-Ray Film Viewer - LED



Unique Features

- Incorporates latest 'Side lit' laser technology for most consistent luminous density
- No drop in light intensity
- Critical components from Japan and Germany ensure reliable performance
- Film insert activated switch
- Continuous brightness control
- Low power consumption
- Sleek design
- Lightest in weight

Model	Viewing Area cm	Luminous Density approx. Lux*	Light Source & Temp	Power Consumption Watts	Power Input	Outer Dim. L x H x B (cm)
Single Film HLX 1	44 x 36	7000	LED 6500 K	15	220-240 VAC	44 x 52 x 2.5
Two Films HLX 2	44 x 72	7000	LED 6500 K	30	220-240 VAC	82 x 52 x 2.5
Three Films HLX 3	44 x 108	7000	LED 6500 K	45	220-240 VAC	120 x 52 x 2.5
Four Films HLX 4	44 x 144	7000	LED 6500 K	60	220-240 VAC	160 x 52 x 2.5

*Luminous density > 10,000 Lux available on request
Product specifications and features are subject to change



Sterling Imaging Solutions
Mumbai 400013 India
E: sterling@sterlingimaging.com | W: www.sterlingimaging.com